



## ELECTROSURGICAL CABLES




**190-230 - Silicone Monopolar Cable 10' (3.0m)**  
Instrument Connector: Screw Type  
Generator Connector: Standar Bovie® 8 mm Plug  
Compatibility: Valleylab, Conmed,Bovie




**190-240 - Silicone Monopolar Endoscopic Cable 10' (3.0m)**  
Instrument Connector: 4 mm Female Fitting  
Generator Connector: Standar Bovie® 8 mm Plug  
Compatibility: Valleylab, Conmed,Bovie




**190-250 - Silicone Monopolar Endoscopic Cable 10' (3.0m)**  
Instrument Connector: 4 mm Female Fitting  
Generator Connector: Standar 4 mm Banana Plug  
Compatibility: Most ES Generators



**190-260 - Silicone Diathermy Cable 10' (3.0m)**  
Instrument Connector: 4 mm Male Fitting(European Diathermy)  
Generator Connector: Standar 4 mm Banana Plug  
Compatibility: Most ES Generators




**190-270 - Silicone Diathermy Cable 10' (3.0m)**  
Instrument Connector: 4.8 mm Female Fitting(English Diathermy)  
Generator Connector: Standar Bovie® 8 mm Plug  
Compatibility: Valleylab, Conmed,Bovie



**190-280 - Silicone Diathermy Cable 10' (3.0m)**  
Instrument Connector: 4.8 mm Female Fitting(English Diathermy)  
Generator Connector: Standar 4 mm Banana Plug  
Compatibility: Most ES Generators




## ELECTROSURGICAL CABLES




**190-100 - Silicone Bipolar Cable 10' (3.0m)**  
Instrument Connector: USA 2-Pin Plug  
Generator Connector: 4 mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus & Most Bipolar Generator



**190-110 - Silicone European Bipolar Cable 10' (3.0m)**  
Instrument Connector: European Flat Plug  
Generator Connector: 4 mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus & Most Bipolar Generator



**190-120 - Silicone European Bipolar Cable 10' (3.0m)**  
Instrument Connector: European Flat Plug  
Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Berchtold, Wolf



**190-130 - Silicone Bipolar Cable 10' (3.0m)**  
Instrument Connector: USA 2-Pin Plug  
Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Berchtold, Wolf




**190-140 - Silicone Bipolar Cable 10' (3.0m)**  
Instrument Connector: European Flat Plug  
Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens




**190-150 - Silicone Bipolar Cable 10' (3.0m)**  
Instrument Connector: USA 2-Pin Plug  
Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens



**190-160 - Silicone European Bipolar Cable 10' (3.0m)**  
Instrument Connector: European Flat Plug  
Generator Connector: Pin Distance 29 mm

Compatibility: Valleylab, EMC, Dolley



**190-170 - Silicone Bipolar Cable 10' (3.0m)**  
Instrument Connector: USA 2-Pin Plug  
Generator Connector: Pin Distance 29 mm

Compatibility: Valleylab, EMC, Dolley



**190-180 - Silicone Bipolar Cable 10' (3.0m)**  
Instrument Connector: Bipolar Artery Sraier  
Generator Connector: 4 mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus & Most Bipolar Generator