




DISPOSABLE BIPOLAR FORCEPS




Disposable Bipolar Cable
190-100D




Jeweler #5
210-140




Adson
210-200



Gerald (Str)
210-300



Jansen Bayonet Irrigation
210-360IR




Jansen Bayonet
210-370

Forceps are supplied with 3m cable.
Available in Sterile/ Non-sterile packaging.
Sterilization by Ethylene Oxide (EO).



More designs are available in Disposable Bipolar Forceps.
Use 210-series with any regular Bipolar Forceps for
conversion into Disposable Bipolar Forceps.



DISPOSABLE BIPOLAR FORCEPS



Jeweler #5
410-110



Straight Forceps
410-140



Bayonet Forceps
410-180

Forceps are supplied with 3m cable.
Available in Sterile/ Non-sterile packaging.
Sterilization by Ethylene Oxide (EO).

More designs are available in Disposable Bipolar Forceps.
Use 410-series with any regular Bipolar Forceps for
conversion into Disposable Bipolar Forceps.



Straight Forceps (Disposable without cable)
510-140



Straight Ang (Disposable without cable)
510-142