

Haemoband


Rapid
Simple
Efficient

Haemoband
SURGICAL





Rapid. Sim



The multi-ligator tip is pre-loaded with four rubber bands

Disposable and low cost, the Haemoband offers a very attractive alternative to current devices by providing an easy to use, ergonomically designed product, pre-loaded with four bands.

- Low cost
- Solo operation
- Multi-application
- Rapid
- Accurate


A slim 120mm insertion probe

Developed with and by surgeons, the unit has been designed to optimise, simplify and substantially cost-reduce the ligation process.

- Reduced patient discomfort
- Single-use
- Fully disposable unit
- Clinically efficient

ple. Efficient.

The Haemoband Multi-Ligator is a new device that facilitates the easy Rubber Band Ligation (RBL) of haemorrhoids.



An innovative dual function hand grip

To use the multi-ligator, position as required and gently squeeze the trigger to induce suction. Continue to pull the trigger to its final position to apply the band, then release the trigger to reduce suction and free the haemorrhoid. The multi-ligator is automatically reloaded with the next band ready for immediate use.