

LifeGuard™ Safety Infusion Sets


BENEFITS:

- > Low-profile for patient comfort
- > Universal luer locks accommodate all injection caps
- > Latex free
- > DEHP-free
- > Soft wings designed for minimal skin irritation
- > Non-siliconized needle helps prevent dislodgement
- > One-handed clamps for positive pressure technique
- > Needleless access connector on Y-Site models

Safe has never been this easy—raise needle cover ninety degrees; hold base firmly, lift wings straight up.

“LifeGuard™ Safety Infusion Sets help us in our commitment to safe delivery of chemotherapy to our patients. The safety device is easy to activate and its low profile adds to patient comfort.”

—Bart Magowan, RN, BSN, OCN; West Clinic, Memphis, TN



LifeGuard™ NAC Y-Site Models	Quantity	Part #
19 Ga X .75" with needleless access connector Y-Site	20 per box	LG-19-75NY
19 Ga X 1" with needleless access connector Y-Site	20 per box	LG-19-100NY
20 Ga X .50" with needleless access connector Y-Site	20 per box	LG-20-50NY
20 Ga X .75" with needleless access connector Y-Site	20 per box	LG-20-75NY
20 Ga X 1" with needleless access connector Y-Site	20 per box	LG-20-100NY
20 Ga X 1.5" with needleless access connector Y-Site	20 per box	LG-20-150NY
22 Ga X .50" with needleless access connector Y-Site	20 per box	LG-22-50NY
22 Ga X .75" with needleless access connector Y-Site	20 per box	LG-22-75NY
22 Ga X 1" with needleless access connector Y-Site	20 per box	LG-22-100NY

LifeGuard™ Non Y-site Models	Quantity	Part #
19 Ga X .75"	20 per box	LG-19-75
19 Ga X 1"	20 per box	LG-19-100
20 Ga X .50"	20 per box	LG-20-50
20 Ga X .75"	20 per box	LG-20-75
20 Ga X 1"	20 per box	LG-20-100
20 Ga X 1.5"	20 per box	LG-20-150
22 Ga X .50"	20 per box	LG-22-50
22 Ga X .75"	20 per box	LG-22-75
22 Ga X 1"	20 per box	LG-22-100

USA > One Horizon Way, Manchester, GA 31816 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
International > Mechelsesteenweg 251, B-1800 Vilvoorde, Belgium > tel: +32-2 252 12 02 > fax: +32-2 252 14 08

www.angiodynamics.com

©2007 AngioDynamics, Inc. pn 0191M rev: B

