


LifePort® Infusion Sets

BENEFITS:

- > Easy to grip, place and control
- > Easy to use with all implantable ports
- > Larger grip handle for secure control
- > Lipid-compatible
- > Latex-free
- > Height-adjustable wings improve needle support and facilitate securing the device to the patient regardless of the patient's size
- > Easy to secure
- > Anti-coring needle
- > One-handed clamp

One Device for All Ports

RITA's LifePort® Infusion Sets are designed for use with all implantable ports. Designed for clinicians, LifePort® Infusion Sets offer large grip wings for secure and maximum control.



LifePort® Infusion Sets With Y-Site/Quantity	Quantity	Part #
19 Ga x 1"	20 per box	LPS 3006
20 Ga x 1.5"	20 per box	LPS 3017
20 Ga x 1"	20 per box	LPS 3009
20 Ga x .75"	20 per box	LPS 3023
22 Ga x 1"	20 per box	LPS 3007
22 Ga x .75"	20 per box	LPS 3025

LifePort® Infusion Sets Without Y-Site/Quantity	Quantity	Part #
19 Ga x 1"	20 per box	LPS 3004
20 Ga x .75"	20 per box	LPS 3010
20 Ga x 1"	20 per box	LPS 3008
22 Ga x 1"	20 per box	LPS 3005

RITA 90° Non-Coring Needles	Quantity	Part #
20 Ga x .75"	10 per box	36860
20 Ga x 1.25"	10 per box	36861

USA > One Horizon Way, Manchester, GA 31816 > tel: 800.472.5221 or 706.846.3126 > fax: 877.467.9501 or 706.846.3146
 International > Mechelsesteenweg 251, B-1800 Vilvoorde, Belgium > tel: +32-2 252 12 02 > fax: +32-2 252 14 08